	Vorkörner Inical bur	Allportbohrer Allport-bur	Lagerfräser Ablative burs		Spezialbohrer Special bur	
Fig.	186RF	I41RF	TC21X	TC21X	BL203	BL203
Schaft Shank <sup>1</sup>	204	204	205	205	206	206
Größe Size <sup>2</sup>	018	031	040	060	045	070
Länge Length mm	12,0	0,3	14,0	14,0	14,0	14,0

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
<sup>1</sup> 204=RA, 206=RA XL    <sup>2</sup> Largest working part diameter in 1/10 mm

## Balloon-Lift-Control

Balloon-Assisted Sinus Floor Elevation

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
- Minimally invasive lifting of internal sinus
- Optimal use of existing bone substance independent of vertical height
- Safe, defined and simple access to the maxillary sinus mucosa
- Complication free lifting by up to approx. 10 mm through extended balloon
- The procedure is easy to pick up and handle
- The use of the procedure saves time and costs
- Suitable to prepare the jaw bone for the implantation of all current systems

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
The balloon catheter can be ordered separately under the following order number 39BLC0K.

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Balloon-Lift-Control (BLC) is a very safe procedure used for the pre-augmentative stripping of the maxillary sinus mucosa. It is minimally invasive and extremely easy to pick up. Independently of the vertical height of the osseous sinus floor, the mucosa can be lifted by more than 10 mm. BLC consists of two osteotomes with the help of which a balloon catheter is being inserted into the maxillary sinus after impression of the weakened floor. Subsequently the balloon is extended to the required size using radio-opaque liquid. The necessary volume equals the volume of augmentation required. This way, sinus floor augmentation becomes a procedure that can be carried out easily and quickly in all dental practices that specialise in implants without leading to the typical complications associated with conventional sinus lifts.

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During the use of the BLC system, this step-by-step approach should be followed:

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The use of the initial bur combined with the respective shell prevents the special bur from slipping. If required, the allport bur or one of the ablative burs can refresh and smooth the bone.

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
Depending on the existing bone and the size of the implant to be inserted at a later date, now is the time to choose the guidance shell, the special bur and the osteotome with the larger (7 mm) or smaller (4.5 mm) diameter. Now the appropriate guidance shell is fastened onto the guidance tool, which is then held in front of the prepared alveolar crest.

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Subsequently and in order to thin the sinus floor you use the special bur to drill through the guidance shell into the alveolar crest until the thickness of the remaining bone measures 1 mm. The perfectly matching tools ensure that this pre-determined height of the remaining bone is maintained.

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Afterwards, the pointed tip of the osteotome is evenly positioned in the pre-shaped canal and the mandrin is inserted.

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
By slightly tapping on the mandrin grip, the compacta on the maxillary sinus side is fractured in a controlled fashion. Prior to this, the part of the mandrin that is now rising above the top of the shaft was accurately adjusted to the height of the remaining bone.

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Removal of the mandrin from the osteotome. Prior to inserting the balloon catheter the air must be evacuated i.e. it must be pre-filled with fluid (see the detailed description of the balloon catheter), then it is inserted into the guiding canal of the osteotome as far as it will go. In doing so, the blue extension at the end of the catheter shaft is pushed right down into the deepening of the osteotome grip.

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This is followed by a repeated (carried out at least 5 times) slow and controlled extension of the balloon with a maximum of approximately 3 ml of fluid (Ultravist® 240 or physiological, 0.9 % NaCl solution). This procedure ensures the stripping and lifting of the maxillary sinus mucosa. If required, a respective X-ray image can be generated with the extended balloon positioned in situ.

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At the end, the balloon is emptied again and removed with the osteotome. Then a spoon is used to insert the exact volume of the augmentation, which before served as the fluid to extend the balloon, step by step through the drilled hole into the newly established sinus area.