
Crestal-Lift-Control

CE 0044

Fig.	186RF	859	TC21X	CL019	CL020	CL001
Shank ¹	204	204	205	206	206	206
Size ²	018	018	040	020	030	020
Length mm	12.0	10.0	14.0	2.5	2.5	13.0

¹ 204=RA, 205=RA L, 206=RA XL ² Largest working part diameter in 1/10 mm

Crestal-Bohrer

Crestal drill
Fresa crestal


Fig.	CL002	CL003	CL004	CL005	CL006	CL007
Shank ¹	206	206	206	206	206	206
Size ²	028	031	033	036	038	041
Length mm	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0

¹ 206=RA XL ² Largest working part diameter in 1/10 mm

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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This system for the internal sinus lift facilitates a simple, safe augmentation of the sinus floor. Elevation occurs during the drilling process. The stop sleeve system is coordinated with special instruments to prevent the membrane from being injured or punctured. In addition to the especially atraumatic design of the Crestal drill with its 4 cutting edges and the concave head for safely forming a conical bone flap, this Crestal drill is also ideal for collecting bone chips.

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
First make an analysis of the residual bone using an x-ray image. After removing the mucous membrane, you use the pilot drill and a suitable stop sleeve for the initial drilling, which ends approx. 2 mm before the sinus. When the thickness of the remaining bone is about 6 mm, you use stop sleeve CL016 to drill 4 mm into the bone and leave the desired 2 mm before the maxillary sinus. During the entire procedure, the area must be cooled with a saline solution. The pilot bur is operated at max. 1,000-1,500 rpm.

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
Now increase the diameters used in the Crestal drill (up to the diameter of the implant, see manufacturer information) using the 6 mm stop sleeve CL014 for the final drilling. This final drilling opens the access to the sinus. The atraumatic design of the Crestal drill with its 4 cutting edges and the concave head allows a conical bone flap to be made. This automatically lifts the membrane by approx. 2 mm. The area must be cooled with a saline solution during the entire procedure and in addition, you should work intermittently. The Crestal drill should be operated at max. 800 rpm. The Crestal drill is ideal for collecting bone chips. Operate at max. 50 rpm.

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
Use the depth gauge to test the mobility of the bone flap. In connection with the appropriate stop sleeve (here CL013, for example), the membrane can be lifted by max. one more millimeter with the depth gauge if needed. If the membrane needs to be lifted by more than 3 mm, you can use the balloon expander (Art. No.: BLC0K) from the Bone Management® System Balloon Lift Control (Art. No.: BBA00) to continue to separate and lift the sinus membrane slowly and under controlled conditions.

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
Then use the spoon to gradually fill the bone material and/or bone chips through the drill hole into the newly created sinus cavity.

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Use the bonegraft condenser to distribute the augmentation material evenly in the sinus cavity that has been created. More bone graft material can then be added if needed. The bone graft spreader should be operated at max. 10 rpm.

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- Safe and fast sinus floor elevation during the drilling process
- Coordinated stop sleeve system prevents injuring and puncturing the membrane
- Atraumatic design of the Crestal drill and depth gauge
- Bone graft spreader optimizes the distribution and compression of the augmentation material

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General instructions:

- All products delivered are unsterile, therefore, before initial and each further potential application, products should be treated (cleaning/disinfection/sterilization).

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